

NanoScience

Tubney Woods, Abingdon, OX13 5QX UK

Tel: +44 (0) 1865 393200

Email: nanoscience@oxinst.com

nanoscience.oxinst.com



59-A60393 50mm LOS Port

LINE No.	RT	RT-70K	70K (PT1)	ATTENUATOR	70K-4K	4K (PT2)	ATTENUATOR	4K-STILL	STILL	ATTENUATOR	STILL-C/PLT	C/PLATE	ATTENUATOR	C/PLT-M/C	M/CHAMBER
1	SMA	UT085 SS/BeCu	SMA	10 dB	UT085 SS/BeCu	SMA	20 dB	UT085 NbTi/NbTi	SMA	06 dB	UT085 NbTi/NbTi	SMA	-----	UT085 NbTi/NbTi	SMA
2	SMA	UT085 SS/BeCu	SMA	10 dB	UT085 SS/BeCu	SMA	20 dB	UT085 NbTi/NbTi	SMA	06 dB	UT085 NbTi/NbTi	SMA	-----	UT085 NbTi/NbTi	SMA
3	SMA	UT085 SS/BeCu	SMA	10 dB	UT085 SS/BeCu	SMA	20 dB	UT085 NbTi/NbTi	SMA	06 dB	UT085 NbTi/NbTi	SMA	-----	UT085 NbTi/NbTi	SMA
4	SMA	UT085 SS/BeCu	SMA	10 dB	UT085 SS/BeCu	SMA	20 dB	UT085 NbTi/NbTi	SMA	06 dB	UT085 NbTi/NbTi	SMA	-----	UT085 NbTi/NbTi	SMA
5	SMA	UT085 SS/BeCu	SMA	10 dB	UT085 SS/BeCu	SMA	20 dB	UT085 NbTi/NbTi	SMA	06 dB	UT085 NbTi/NbTi	SMA	-----	UT085 NbTi/NbTi	SMA
6	SMA	UT085 SS/BeCu	SMA	10 dB	UT085 SS/BeCu	SMA	20 dB	UT085 NbTi/NbTi	SMA	06 dB	UT085 NbTi/NbTi	SMA	-----	UT085 NbTi/NbTi	SMA
7	SMA	UT085 SS/BeCu	SMA	10 dB	UT085 SS/BeCu	SMA	20 dB	UT085 NbTi/NbTi	SMA	06 dB	UT085 NbTi/NbTi	SMA	-----	UT085 NbTi/NbTi	SMA
8	SMA	UT085 SS/BeCu	SMA	10 dB	UT085 SS/BeCu	SMA	20 dB	UT085 NbTi/NbTi	SMA	06 dB	UT085 NbTi/NbTi	SMA	-----	UT085 NbTi/NbTi	SMA

Please note: Where no attenuator is fitted the coax is longer and connects between both bulkheads.

59-A60398 50mm LOS Port

LINE No.	RT	RT-70K	70K (PT1)	ATTENUATOR	70K-4K	4K (PT2)	ATTENUATOR	4K-STILL	STILL	ATTENUATOR	STILL-C/PLT	C/PLATE	ATTENUATOR	C/PLT-M/C	M/CHAMBER
1	SMA	UT085 CuNi/CuNi	SMA	10 dB	UT085 CuNi/CuNi	SMA	20 dB	UT085 NbTi/NbTi	SMA	06 dB	UT085 NbTi/NbTi	SMA	-----	UT085 NbTi/NbTi	SMA
2	SMA	UT085 CuNi/CuNi	SMA	10 dB	UT085 CuNi/CuNi	SMA	20 dB	UT085 NbTi/NbTi	SMA	06 dB	UT085 NbTi/NbTi	SMA	-----	UT085 NbTi/NbTi	SMA
3	SMA	UT085 CuNi/CuNi	SMA	10 dB	UT085 CuNi/CuNi	SMA	20 dB	UT085 NbTi/NbTi	SMA	06 dB	UT085 NbTi/NbTi	SMA	-----	UT085 NbTi/NbTi	SMA
4	SMA	-----	SMA	---	-----	SMA	---	-----	SMA	---	-----	SMA	-----	-----	SMA
5	SMA	UT085 CuNi/CuNi	SMA	10 dB	UT085 CuNi/CuNi	SMA	20 dB	UT085 NbTi/NbTi	SMA	06 dB	UT085 NbTi/NbTi	SMA	-----	UT085 NbTi/NbTi	SMA
6	SMA	UT085 CuNi/CuNi	SMA	10 dB	UT085 CuNi/CuNi	SMA	20 dB	UT085 NbTi/NbTi	SMA	06 dB	UT085 NbTi/NbTi	SMA	-----	UT085 NbTi/NbTi	SMA
7	SMA	UT085 CuNi/CuNi	SMA	10 dB	UT085 CuNi/CuNi	SMA	20 dB	UT085 NbTi/NbTi	SMA	06 dB	UT085 NbTi/NbTi	SMA	-----	UT085 NbTi/NbTi	SMA
8	SMA	-----	SMA	---	-----	SMA	---	-----	SMA	---	-----	SMA	-----	-----	SMA

Please note: Where no attenuator is fitted the coax is longer and connects between both bulkheads. Line 4 and 8 no coax fitted, bulkhead SMA only.

59-A212343 50mm LOS Port

LINE No.	RT	RT-70K	70K (PT1)	ATTENUATOR	70K-4K	4K (PT2)	ATTENUATOR	4K-STILL	STILL	ATTENUATOR	STILL-C/PLT	C/PLATE	ATTENUATOR	C/PLT-M/C	M/CHAMBER
1	SMA	UT085 SS/SS	SMA	10 dB	UT085 SS/SS	SMA	20 dB	UT085 SS/SS	SMA	06 dB	UT085 SS/SS	SMA	-----	UT085 SS/SS	SMA
2	SMA	UT085 SS/SS	SMA	10 dB	UT085 SS/SS	SMA	20 dB	UT085 SS/SS	SMA	06 dB	UT085 SS/SS	SMA	-----	UT085 SS/SS	SMA
3	SMA	UT085 SS/SS	SMA	10 dB	UT085 SS/SS	SMA	20 dB	UT085 SS/SS	SMA	06 dB	UT085 SS/SS	SMA	-----	UT085 SS/SS	SMA
4	SMA	UT085 SS/SS	SMA	10 dB	UT085 SS/SS	SMA	20 dB	UT085 SS/SS	SMA	06 dB	UT085 SS/SS	SMA	-----	UT085 SS/SS	SMA
5	SMA	UT085 SS/SS	SMA	10 dB	UT085 SS/SS	SMA	20 dB	UT085 SS/SS	SMA	06 dB	UT085 SS/SS	SMA	-----	UT085 SS/SS	SMA
6	SMA	UT085 SS/SS	SMA	10 dB	UT085 SS/SS	SMA	20 dB	UT085 SS/SS	SMA	06 dB	UT085 SS/SS	SMA	-----	UT085 SS/SS	SMA
7	SMA	UT085 SS/SS	SMA	10 dB	UT085 SS/SS	SMA	20 dB	UT085 SS/SS	SMA	06 dB	UT085 SS/SS	SMA	-----	UT085 SS/SS	SMA
8	SMA	UT085 SS/SS	SMA	10 dB	UT085 SS/SS	SMA	20 dB	UT085 SS/SS	SMA	06 dB	UT085 SS/SS	SMA	-----	UT085 SS/SS	SMA

Please note: Where no attenuator is fitted the coax is longer and connects between both bulkheads.

Please do not share this document with anyone outside of your organisation

NanoScience

Tubney Woods, Abingdon, OX13 5QX UK

Tel: +44 (0) 1865 393200

Email: nanoscience@oxinst.com

nanoscience.oxinst.com

**59-A212451 50mm LOS Port**

LINE No.	RT	RT-70K	70K (PT1)	ATTENUATOR	70K-4K	4K (PT2)	ATTENUATOR	4K-STILL	STILL	ATTENUATOR	STILL-C/PLT	C/PLATE	ATTENUATOR	C/PLT-M/C	M/CHAMBER
1	SMA	UT085 SS/SS	SMA	10 dB	UT085 SS/SS	SMA	20 dB	UT085 SS/SS	SMA	06 dB	UT085 SS/SS	SMA	-----	UT085 SS/SS	SMA
2	SMA	-----	SMA	---	-----	SMA	---	-----	SMA	---	-----	SMA	-----	-----	SMA
3	SMA	UT085 SS/SS	SMA	10 dB	UT085 SS/SS	SMA	20 dB	UT085 SS/SS	SMA	06 dB	UT085 SS/SS	SMA	-----	UT085 SS/SS	SMA
4	SMA	-----	SMA	---	-----	SMA	---	-----	SMA	---	-----	SMA	-----	-----	SMA
5	SMA	UT085 SS/SS	SMA	10 dB	UT085 SS/SS	SMA	20 dB	UT085 SS/SS	SMA	06 dB	UT085 SS/SS	SMA	-----	UT085 SS/SS	SMA
6	SMA	-----	SMA	---	-----	SMA	---	-----	SMA	---	-----	SMA	-----	-----	SMA
7	SMA	UT085 SS/SS	SMA	10 dB	UT085 SS/SS	SMA	20 dB	UT085 SS/SS	SMA	06 dB	UT085 SS/SS	SMA	-----	UT085 SS/SS	SMA
8	SMA	-----	SMA	---	-----	SMA	---	-----	SMA	---	-----	SMA	-----	-----	SMA

Please note: Where no attenuator is fitted the coax is longer and connects between both bulkheads. Line 2, 4, 6 and 8 no coax fitted, bulkhead SMA only.

Please do not share this document with anyone outside of your organisation